



# FilterTEF 0.2 micron

Validation Guide



SISTEMI DI FILTRAZIONE

Via Madre Teresa, 22 - 20851 – Lissone (MB), Italy

Tel. +39.039.2145244 - Fax. +39.039.2145245

[www.ipmfilters.com](http://www.ipmfilters.com)

## ➤ Summary

- General information filter cartridges FilterTEF	pag. 3
- Technical characteristics FilterTEF rated 0.2 micron	pag. 4
- Bacteria Removal efficiency	pag. 5
- Challenge Test filter cartridge type TEF J05 002 SI-PH	pag. 6/ 7
- Correlation diffusion vs retention cartridge type TEF 710 002 SI -PH	pag. 8/ 9
- Correlation diffusion vs retention cartridge type TEF 720 002 SI – PH	pag. 10/11
- Biological and extractables tests	pag. 12
- Bacterial endotoxins	pag. 13
- Chemical compatibility Extractables as per USP 661	pag. 14
- Biological safety test Test “in vitro” class V USP 87	pag. 15
Test “in vitro” classe VI USP 88	pag. 16
- Physical Proprieties Hydraulic stress test	pag. 17
Thermal stress test	pag. 18
- Delta P Vs water flow rate	pag. 19
- Attached: Certificate of conformity	pag. 20

## ➤ General information filter cartridges FilterTEF

This document give the necessary information for the use of *FilterTEF* rated 0,2 micron in pharmaceutical applications. .

FilterTEF filter cartridges 0,2 micron rated, are manufactures as per the procedures used for production of sterilizing grade filters.

Materials used for the assembly of the cartridges, both the teflon (PTFE) membrane and the Polypropylene used for the hardware, had been selected for the specific application meeting the current rules.

The assembly is made by thermo-bonding without any adhesive.

*FilterTEF* filter cartridges rated 0,2 micron are manufactured into a clean room with controlled contamination and tested prior to the shipment by the internal IPM labs in order to meet the high standards and the quality assurance as per internal procedures.

The identification of the FilterTEF filter cartridges is made per the following:

TEF	7	10	002	S	I	PH
	CARTRIDGE TYPE	LENGTH	FILTRATION RATING	GASKETS O-RINGS	INSERT	VERSION
Product identification: <b>TEF</b>	D = DOE (double open end) 2 = bayonet 2.226/flat 3 = 2.222/flat 5 = 3 bayonet 2.222/fin 7 = bayonet 2.226/fin 8 = 2.222/fin	10 = 10" 20 = 20" 30 = 30" 40 = 40"	001 = 0.1 µm 002 = 0.2 µm 004 = 0.45 µm	B = Buna N E = EPR S = Silicone V = Viton X = other  VITON® is a registered trade name of E. I. duPont de Nemours & Co. Inc	I = with SS insert (*)  0 = no SS insert  (*) When filter has to be sterilized by steam	PH= Pharma version
	J = Junior Cartridge	02 = 2,5" 05 = 5"				

Should you need any more info please contact IPM technical support tel. +39 039 21 45 244 or visit the site [www.ipmfilters.com](http://www.ipmfilters.com)

➤ **Technical characteristics FilterTEF rated 0.2 micron**

<b>Cartridges Materials:</b>				
<i>Membrane</i>	Teflon (PTFE)			
<i>Support Layers</i>	Polypropylene (PP)			
<i>Cartridges supprts</i>	Polypropylene (PP)			
<i>O-rings</i>	Silicon			
<b>Technical Characteristics of the cartridge:</b>				
<i>Filtration Area</i>	0.78 m <sup>2</sup> per 10" modules			
<i>Pore size</i>	0.2 µm			
<i>Length</i>	Available cartridges as per the following dimensions: 5" JUNIOR (125 mm), 10" (250 mm), 20" (500 mm), 30" (750 mm)			
<b>Retention characteristics of the cartridges:</b>				
<i>Bacteria Retention</i>	Retention of 10 <sup>7</sup> CFU/cm <sup>2</sup> <i>Brevundimonas diminuta</i> (ATCC 19146), as per ASTM F838-05.			
<b>Mechanic Characteristics</b>				
<i>Resistance at Hydraulic Stress</i>	<i>FilterTEF 0,2 µm</i> filter cartridges shown no diffusion and physical changes when stressed by:			
	Water Pressure (bar)	N° strokes	Water temperature (°C)	Direction of the flow
	5,0	50	20±5	Forward
	2,0	50	20±5	Reverse
<i>Resistance at Thermal Stress</i>	<i>FilterTEF 0,2 µm</i> filter cartridges shown no diffusion and physical changes when stressed by:			
	Temperature (°C)	Test	Pressure Test (bar)	Cycle <sup>1</sup> time (hours)
	121		1.15	1
				N° of cycles (in autoclave)
				130
<b>Extractables</b>				
<i>Extractables</i>	<i>FilterTEF 0,2 µm</i> filter cartridges meet USP<661>.			
<i>Conformity of materials to the rules for food contact</i>	Materials of the <i>FilterTEF 0,2 µm</i> cartridge meet the requirements for food contact as per Reg. N° 1935/2004/EEC and "Code of Federal Regulations" FDA, 21 # 171-199.			
<b>Biological Reactivity</b>				
<i>Material toxicity</i>	All the materials of the <i>FilterTEF 0,2 µm</i> cartridges meet the rules as per USP<87> ("Biological Reactivity Test in vitro for Class VI Plastics") e USP<88> ("Biological Reactivity Test in vivo for Class VI Plastics – Safety Test").			
<b>Release of particles and/or Endotoxins</b>				
<i>Bacteria endotoxins</i>	< 0,125 EU/ml of aqueous extract as per USP<85>.			
<i>No fibers releasing</i>	<i>FilterTEF 0,2 µm</i> cartridges satisfy the rules as per USP<788>			

<sup>1</sup> Total time considering a warming up rising of 15 min, 30 min at the test temperature and of 15 min. for room temperature cooling down

## ► Bacteria Removal efficiency

Several Authorities as FDA and EMA fix the standards for show that the sterilizing action of a filter has to be evidenced by tests before to be used in the pharmaceutical application. A filter is cannot be considered sterilizing grade only because the dimension of its pores.

In order to consider a filter as sterilizing grade it must challenged with an upstream concentration not less of  $10^7$  CFU/cm<sup>2</sup> of the bacterial *Brevundimonas diminuta*. It must get a sterile downstream.

Correlation between the microbiological retention and a non destructive integrity test is an other very important aspect of the validation of a filter cartridge pharmaceutical grade.

The Official Authorities had set up a procedure for control the filter both prior to the use and after the use.

On the followings you can find the Challenge tests results to confirm that the membrane used for the filter cartridges *FilterTEF 0,2 μm* can have full retention a bacteria charge of  $10^7$  CFU/cm<sup>2</sup> di *Brevundimonas diminuta* according to the ASTM F838-05 methodology.

You can also find the correlation between the retention and the diffusion test.

This is important to set the maximum diffusion flow acceptance limit to warranty the filter will get a sterile effluent.

On the followings tables are reported the correlation between the diffusion test (non destructive test) and the bacteria challenge (destructive test) for filter cartridges *FilterTEF 0,2 μm* Part Numbered: TEF J05 002 SI - PH ; TEF 710 002 SI - PH ; TEF 720 002 SI - PH.

### ➤ Challenge Test - Filter cartridge TEF J05 002 SI - PH da 0.2 micron

20 (twenty) cartridges from the standard production had been tested. Any of them had been tested first to an integrity test by diffusion and then bacteria challenged by *Brevundimonas diminuta*  $>10^7$  CFU/cm<sup>2</sup>.


Table 1 shows the Tests results also shown by graphic on Fig 1.

As per the results the filter cartridge TEF J05 002 SI – PH - 0.2 micron is conform to the acceptance requirements.

Lot	Diffusion at 690 mbar (ml/min) **	Downstream Sterile	Titre Reduction
6818/22	4,0	Si	$> 1,6 \times 10^8$
6818/25	4,8	Si	$> 1,5 \times 10^8$
6117/06	5,1	Si	$> 1,6 \times 10^8$
6118/05	5,3	Si	$> 1,5 \times 10^8$
6117/18	5,4	Si	$> 1,7 \times 10^8$
6118/23	6,0	Si	$> 1,8 \times 10^8$
6777/11	6,5	Si	$> 1,4 \times 10^8$
6771/16	7,0	Si	$> 1,5 \times 10^8$
6512/01	7,3	Si	$> 1,5 \times 10^8$
5943/09	7,3	Si	$> 1,6 \times 10^8$
5836/14	7,4	Si	$> 1,7 \times 10^8$
6132/12	7,5	Si	$> 1,8 \times 10^8$
5981/17	7,6	Si	$> 1,9 \times 10^8$
6091/19	7,6	Si	$> 1,7 \times 10^8$
5643/02	7,9	Si	$> 1,6 \times 10^8$
6843/13	8,0	Si	$> 1,5 \times 10^8$
5643/05	8,0	Si	$> 1,7 \times 10^8$
5643/11	8,0	Si	$> 1,4 \times 10^9$
6777/15	8,5	No	$1,3 \times 10^6$
6777/24	9,0	No	$1,7 \times 10^6$
6818/03	10,5	Si	$> 1,7 \times 10^8$

**Tab. 1. Bacteria challenge results for FilterTEF J05 002 SI - PH - 0.2 µm**

\*\* Tested in 60% IPA / 40% DI Water



**Fig. 1. Correlation between diffusion test and bacteria challenge ASTM F838-05 for filter cartridge FilterTEF J05 002 SI - PH - 0.2 micron**

**Conclusions:**

Based upon the results the filter cartridge TEF J05 002 SI - PH - 0.2 micron meets the acceptance and has to be considered a sterilizing grade filter up to a maximum diffusion value of 8.0 mL/min.

### ► Correlation diffusion Vs cartridge retention TEF 710 002 SI – PH - 0.2 micron

20 (twenty) cartridges from the standard production had been tested. Any of them had been tested first to an integrity test by diffusion and then bacteria challenged by *Brevundimonas diminuta*  $>10^7$  CFU/cm<sup>2</sup>.

Table 2 shows the Tests results also shown by graphic on Fig 2.


As per the results the filter cartridge TEF 710 002 SI – PH - 0.2 micron is conform to the acceptance requirements.

Lot serial N#	Diffusion at 690 mbar ml/min **	Downstream Sterile	Titre Reduction
7095/10	12,0	Si	$> 1,4 \times 10^8$
7095/15	12,3	Si	$> 1,4 \times 10^8$
6827/34	13,0	Si	$> 1,2 \times 10^8$
7095/43	13,3	Si	$> 1,6 \times 10^8$
6869/05	13,4	Si	$> 1,7 \times 10^8$
6726/18	13,7	Si	$> 1,8 \times 10^8$
6726/06	14,0	Si	$> 1,8 \times 10^8$
6366/30	14,1	Si	$> 1,4 \times 10^8$
7095/03	14,3	Si	$> 1,5 \times 10^8$
6566/42	14,3	Si	$> 1,5 \times 10^8$
6827/05	14,4	Si	$> 1,3 \times 10^8$
6726/15	14,5	Si	$> 1,4 \times 10^8$
6366/14	14,6	Si	$> 1,2 \times 10^8$
6566/19	14,6	Si	$> 1,1 \times 10^8$
6366/19	14,9	Si	$> 1,4 \times 10^8$
6827/21	15,0	Si	$> 1,7 \times 10^8$
6566/02	15,0	Si	$> 1,8 \times 10^8$
6366/27	15,0	Si	$> 1,3 \times 10^9$
6827/16	15,5	No	$1,6 \times 10^6$
6869/12	16,0	No	$1,5 \times 10^6$
6366/21	17,5	Si	$> 1,7 \times 10^8$

**Tab. 2. Bacteria challenge results for FilterTEF 710 002 SI - PH - 0.2 µm**

\*\* Tested in 60% IPA / 40% DI Water





**Fig. 2. Correlation between diffusion test and bacteria challenge ASTM F838-05 for filter cartridge FilterTEF 710 002 SI - PH - 0.2 micron**

**Conclusions:**

Based upon the results the filter cartridge TEF 710 002 SI - PH - 0.2 micron meets the acceptance and has to be considered a sterilizing grade filter up to a maximum diffusion value of 15.0 mL/min.

➤ **Correlation diffusion Vs cartridge retention TEF 720 002 SI – PH - 0.2 micron**

20 (twenty) cartridges from the standard production had been tested. Any of them had been tested first to an integrity test by diffusion and then bacteria challenged by *Brevundimonas diminuta*  $>10^7$  CFU/cm<sup>2</sup>.


Table 3 shows the Tests results also shown by graphic on Fig 3.

As per the results the filter cartridge TEF 720 002 SI – PH - 0.2 micron is conform to the acceptance requirements.

Lot	Diffusion at 690 mbar ml/min **	Downstream Sterile	Titre Reduction
5823/14	26,0	Si	$> 1,3 \times 10^8$
5823/23	26,2	Si	$> 1,3 \times 10^8$
4713/23	27,0	Si	$> 1,4 \times 10^8$
4713/12	27,1	Si	$> 1,4 \times 10^8$
5901/05	27,4	Si	$> 1,7 \times 10^8$
5901/16	27,7	Si	$> 1,3 \times 10^8$
4961/12	28,2	Si	$> 1,6 \times 10^8$
5823/24	28,3	Si	$> 1,8 \times 10^8$
5901/22	28,5	Si	$> 1,6 \times 10^8$
6004/11	28,7	Si	$> 1,5 \times 10^8$
4239/19	28,7	Si	$> 1,4 \times 10^8$
5097/03	29,1	Si	$> 1,5 \times 10^8$
4713/18	29,4	Si	$> 1,8 \times 10^8$
4239/15	29,5	Si	$> 1,3 \times 10^8$
4713/19	29,8	Si	$> 1,6 \times 10^8$
5097/32	29,9	Si	$> 1,4 \times 10^8$
5097/17	30,0	Si	$> 1,5 \times 10^8$
6004/15	30,0	Si	$> 1,6 \times 10^9$
4713/32	30,4	No	$1,7 \times 10^6$
5901/30	31,0	No	$1,6 \times 10^4$
5823/04	32,7	Si	$> 1,6 \times 10^8$

**Tab. 3. Bacteria challenge results for FilterTEF 720 002 SI - PH - 0.2 µm**

\*\* Tested in 60% IPA / 40% DI Water



**Fig. 3. Correlation between diffusion test and bacteria challenge ASTM F838-05 for filter cartridge FilterTEF 720 002 SI - PH - 0.2 micron**

**Conclusions:**

Based upon the results the filter cartridge TEF 720 002 SI - PH - 0.2 micron meets the acceptance and has to be considered a sterilizing grade filter up to a maximum diffusion value of 30.0 mL/min.

## Biological and extractables tests

### ➤ Bacteria endotoxins

Test had been run according the Current USP<85> using the Chromogenic Method that gives the endotoxyn amount released by a filter cartridge type *FilterTEF 0,2 micron*.

The sample cartridge had been filled with a-pyrogenic water and ultrasound treated in order to better extract the endotoxins.

The lab test shows no Bacteria endotoxins release by the filter cartridge FilterTEF.

LOT	METHOD	RESULT
5795/03	cUSP<85>	<0.25 EU/mL
5901/17	cUSP<85>	<0.25 EU/mL
4713/29	cUSP<85>	<0.25 EU/mL

*Tabella 4. LAL test Results*

➤ Chemical compatibility and extractables

On the following table the chemical compatibility of *FilterTEF 0.2 micron* challenged with several chemicals. However, in case of concerns, before using *FilterTEF 0.2 micron* with no experienced chemicals, please contact IPM.

<b>CHEMICAL CLASSIFICATION</b>	<b>COMPOUND</b> <i>(Condizioni di prova: 7 giorni a 20°C)</i>	<b>FILTER TEF</b>
ACIDS	Acetic acid, 90%	R
	Acetic acid, 30%	R
	Acetic acid, 10%	R
	Hydrochloric acid, conc. (35%)	R
	Hydrochloric acid, 6N (20%)	R
	Hydrochloric acid, 1N (3.3%)	R
	Nitric Acid, conc. (67%)	R
	Nitric acid, 6N (27%)	R
	Sulfuric acid, conc. (96%)	R
Sulfuric acid, 6N (16%)	R	
ALCOHOLS	Amyl alcohol	R
	Benzyl alcohol	R
	Butanol	R
	Ethanol	R
	Isopropanol	R
	Methanol	R
AROMATIC HYDROCARBONS	Benzene	R
	Toluene	R
	Xylene	R
BASES	Ammonium hydroxyde, 3N (5.7%)	R
	Ammonium hydroxyde, 6N (11.4%)	R
	Potassium hydroxyde, 3N (15%)	R
	Sodium hydroxyde, 3N (11%)	R
	Sodium hydroxyde, 6N (22%)	R
ESTERS	Amyl acetate	R
	Butyl acetate	R
	Cellosolve acetate	R
	Ethyl acetate	R
	Isopropyl acetate	R
	Methyl acetate	R
ETHERS	Ethyl ether	R
	Tetrahydrofuran (THF)	R
	THF in water (50:50)	R
GLYCOLS	Ethylene glycol	R
	Glycerol	R
	Propylene glycol	R
HALOGENATED HYDROCARBONS	Carbon tetrachloride	R
	Chloroform	R
	Ethylene dichloride	R
	Methylene chloride	R
KETONES	Acetone	R
	Cyclohexanone	R
	Methyl ethyl ketone (MEK)	R
	Methyl isobutyl ketone	R
OILS	Cottonseed	R
	Peanut	R
MISCELLANEOUS	Acetonitrile	R
	Dimethylformamide (DMF)	R
	Formaldehyde (37%)	R
	Hexane	R
	Pyridine	R

**Tabella 5. Chemical compatibility FilterTEF da 0.2 micron (R=Recomended, L=Limited Resistance, N=No Recomend)**

➤ Extractables as per USP<661>

Purpose of the following test is to verify the non volatile residuals amount that can be extract by WFI. Test had been run in static condition at 70°C for 24 h.

After the extraction procedure a sample had been dry evaporated to get the Non Volatile Residual (NVR)

LOT	NVR
5795/07	<5 mg
5901/16	<5 mg
4713/28	<5 mg

**Tab. 6. FilterTEF 0.2 micron NVR by la cUSP<661>**

*Conclusions:*

*FilterTEF da 0.2 micron meet NVR by USP<661>.*

## Biological safety test

### ➤ Test “in vitro” class V USP 87

The purpose of this test is to evaluate the biological toxicity of the filter cartridges FilterTEF da 0,2 micron. The sample cartridges, including the o-rings, are treated by the solutions as per procedures described in the current USP<87> in the indicated rapport weigh /volume.

Test result is indicated as a grade of reactivity where 0 is no reactivity and 4 is high reactivity.

Results in the following table:

LOT	TEST	METHOD	RESULT
5795/08	citotossicità	cUSP<87>	Reactivity grade 0
5901/10	citotossicità	cUSP<87>	Reactivity grade 0
4713/25	citotossicità	cUSP<87>	Reactivity grade 0

*Tabella 7. tests by cUSP<87>*

Conclusion:

No in vitro biological reactivity for the filter cartridges FilterTEF da 0,2 micron as per USP<87>.

➤ "in vivo" Class VI plastic materials Safety Test (cUSP<88>)

Purpose of the test is to evaluate the "in vivo" biological reactivity of the filter cartridges FilterTEF as per USP<88>.

Results had been evaluated as meeting the limit described in USP procedures.

<b>LOT</b>	<b>TEST</b>	<b>RESULT</b>
5795/10	Intracutaneous Test	Conform
	Systemic Injection Test	Conform
	Implantation Test	Conform
5901/11	Intracutaneous Test	Conform
	Systemic Injection Test	Conform
	Implantation Test	Conform
4713/22	Intracutaneous Test	Conform
	Systemic Injection Test	Conform
	Implantation Test	Conform

*Tabella 8. Results by cUSP<88>*

Conclusions:

No biological reactivity "in vivo" for the filter cartridges FilterTEF - 0,2 micron by USP<88>.



## Physical Properties

### ➤ Thermal stress test

Purpose of the test is to verify the resistance of the filter cartridge FilterTEF – at repetitive steam sterilizing cycles.

Test had been run at 121°C, 35 min., per cycle.

Sample is tested by diffusion integrity test and tested each 10 cycles by diffusion.

Tab. 9. shows the schedule.

FASE	TEMPERATURE (°C)	PRESSURE (bar)
STERILIZATION	121±1	1.14±0.04
WARMIG UP	121→96	1.12→0.00

*Tab. 2. Sterilizing cycle program*

Diffusion integrity test run at 690 mbar with acceptance limit of di 15,0 mL/min.

Tab.10. Shows the results.

LOT	DIFFUSION BEFORE (mL/min)	DIFFUSION AFTER (mL/min)
6827/02	13,0	14,0
6566/10	14,0	14,5
7095/26	13,5	13,0
6366/05	13,0	14,0
7095/19	12,5	13,5
6726/16	14,5	13,5
6566/12	13,0	14,5
6869/18	14,5	15,0
6566/35	13,0	13,5

*Tab. 3. Thermic stress FilterTEF*

Conclusion:

The filter has high resistance when challenged at several sterilizing cycles.

➤ **Pressure drop Vs water flow rate**

The following table shows the pressure drop Vs. air flow rate at 20±5°C for the filter TEF per 10" module.

In the typical range of use of the filterTEF cartridges the initial Delta P are less than 0,2 bar (flow rate <1200 L/h per 10").




Fig. 2 – Flow rate – Delta P FilterTEF per 10" module

TEF	7	10	002	S	I	PH
Product identification: <b>TEF</b>	<b>CARTRIDGE TYPE</b>	<b>LENGTH</b>	<b>FILTRATION RATING</b>	<b>GASKETS O-RINGS</b>	<b>INSERT</b>	<b>VERSION</b>
	D = DOE (double open end) 2 = bayonet 2.226/flat 3 = 2.222/flat 5 = 3 bayonet 2.222/fin 7 = bayonet 2.226/fin 8 = 2.222/fin  J = Junior Cartridge	10 = 10" 20 = 20" 30 = 30" 40 = 40"  02 = 2,5" 05 = 5"	001 = 0.1 μm 002 = 0.2 μm 004 = 0.45 μm	B = Buna N E = EPR S = Silicone V = Viton X = other  VITON® is a registered trade name of E.I. duPont de Nemours & Co. Inc	I = with SS insert (*)  0 = no SS insert	PH= Pharma version
					(*) When filter has to be sterilized by steam	

attached

➤ **Certificate of Conformity**

Every cartridge for pharma use comes with a CofC like the following fac-simile:

 <p>Via Madre Teresa, 22 20035 Lissone - Italy Tel. +39 039 2140244 Fax +39 039 2140245 E-mail: ipm@ipmfilters.com www.ipmfilters.com</p>		<p><i>Data di produzione</i> Date of production Date de production Fecha de producción</p> <p>01/01/2013</p>	
<p><b>CERTIFICATO DI CONFORMITÀ</b> <b>CERTIFICATE OF CONFORMITY</b> <b>CERTIFICAT DE CONFORMITE</b> <b>CERTIFICADO DE CONFORMIDAD</b></p>		<p>XXXXXXXXXX</p>	
<p><i>Cliente - Customer</i> Client - Cliente</p>		<p>XXXXXXXXXXXXXXXXXX</p>	
<p><i>TIPO DI FILTRO</i> FILTER TYPE TYPE DE FILTRE TIPO FILTRO</p>		<p>Filter TEF</p>	
<p><i>MOD.</i></p>		<p>TEF710.002-S1</p>	
<p><i>LOT.</i></p>		<p>L/N: XXXX</p>	
<p><i>Controllo visivo e dimensionale:</i> Visual and dimensional check: <i>Contrôle visuel et dimensionnel:</i> Inspección visual y control dimensional:</p>		<p>✓</p>	
<p><i>Risultato - Result - Résultat - Resultado</i></p>		<p>ACCEPTED</p>	
<p><i>Test di integrità (diffusione):</i> Integrity test (diffusion): <i>Test d'intégrité (diffusion):</i> Prueba de integridad (difusión):</p>		<p>✓</p>	
<p><i>Limite di accettabilità:</i> Accepted values: <i>Limite de tolérance:</i> Valor límite aceptable:</p>		<p>As per IPM QC procedures</p>	
<p><small>Il test di integrità come sopra specificato, è stato eseguito su ogni singolo cartuccia di questo tipo e risponde alle procedure interne del controllo qualità. IPM certifies that the filter cartridge IPM rigida described above was conducted in ambient controlled. Materiales de construcción: Las componentes del filtro cumplen las especificaciones de integridad suministradas en el listado adjunto de la United States Pharmacopoeia (USP) para Clase VI de plástico a 121°C. These filters also are made from materials listed for total contact usage per Title 21 of the U.S. Code of Federal Regulations (CFR) parts 370-159. Contactar con IPM para mayores informaciones relativas a los materiales de fabricación. Este producto cumple todos los estándares de calidad IPM. El producto no se entrega esterilizado. IPM declara que este certificado de conformidad con instalaciones y se corresponde con los requisitos de su gestión. Todos los procedimientos utilizados se han llevado a cabo con arreglo a sus procedimientos. Los productos y los materiales están libres de defectos.</small></p>			
<p><b>IL CONTROLLO DI QUALITÀ - QUALITY CONTROL</b> <b>CONTRÔLE DE QUALITÉ - CONTROL DE CALIDAD</b></p>			

Fig. 3. Certificate of conformity released by IPM S.r.l.